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The image shows a virtual chemistry laboratory interface. At the top, there is a navigation bar with tabs for '教学' (Teaching), '信息' (Information), '视角' (View), and '其他功能' (Other Functions). To the right of these tabs, there is a user account section labeled '体验账号' (Experience Account) with a refresh icon and the number '0.0'. Below this, a tooltip indicates '按 [Esc] 即可退出全屏模式' (Press [Esc] to exit full-screen mode). The main area displays a 3D-rendered laboratory setup. In the center, there is a piece of equipment with two gauges and a blue base, connected to a glass flask on a stand. To the right, there is a table with a beaker, a small yellow box, and a graduated cylinder. A large white bottle labeled '蒸馏水' (Distilled Water) and a smaller bottle labeled '稀氨水' (Dilute Ammonia Water) are also visible. The background shows a typical laboratory environment with blue cabinets and a sink. On the right side, there is a circular control panel with icons for '换肤' (Change Skin), '反馈' (Feedback), '退出' (Exit), '帮助' (Help), and '全屏' (Full Screen). A watermark '南京师范大学教师教学发展中心版权所有' (Copyright © Nanjing Normal University Teacher Education Development Center) is visible across the center of the image.

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录制：EV Capture软件

The screenshot displays the EV Capture software interface. At the top, there is a blue header bar with the text "EV录屏" on the left and a "设置" (Settings) button circled in red on the right. Below the header, the interface is divided into several sections:

- 常规 (General):** Includes buttons for "本地录制" (Local Recording) and "在线直播" (Live Streaming).
- 视频 (Video):** Contains radio buttons for "全屏录制" (Full Screen Recording), "选区录制" (Region Recording), "摄像头录制" (Camera Recording), and "不录视频" (Do Not Record Video).
- 音频 (Audio):** Contains radio buttons for "仅麦克风声音" (Microphone Only), "仅系统声音" (System Only), "麦和系统声音" (Microphone and System), and "不录音频" (Do Not Record Audio).
- 辅助工具 (Auxiliary Tools):** Includes buttons for "图片水印" (Image Watermark), "文字水印" (Text Watermark), "嵌入摄像头" (Embed Camera), and "本地直播" (Local Live Streaming).
- 预览窗口 (Preview Window):** Shows a desktop environment with a file explorer window open. A red arrow points to the "设置" button in the top right corner of the preview window.
- 底部控制栏 (Bottom Control Bar):** Features a play button, a stop button, a timer showing "00:00:00", and a volume control icon.

A watermark "南京大学教师教学发展中心版权所有" is visible diagonally across the interface.

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录制：EV Capture软件

The image shows the EV录屏 (EV Capture) software interface. The main window has a blue header with the text "EV录屏" and a user icon on the left, and "设置" (Settings) with window control icons on the right. Below the header, there are two main tabs: "本地录制" (Local Recording) and "在线直播" (Live Streaming). The "本地录制" tab is active. Underneath, there are sections for "视频" (Video) and "音频" (Audio). A "设置" (Settings) dialog box is open in the foreground, showing various configuration options. The settings dialog has a blue header with "设置" and a close button. It is divided into several sections: "通用设置" (General Settings) with a text input for "地址密钥" (Address Key); "本地存盘" (Local Storage) with a "快捷键" (Shortcuts) section containing: "开始/暂停" (Start/Pause) set to Ctrl + F1, "停止" (Stop) set to Ctrl + F2, and "显示主窗口" (Show Main Window) set to Ctrl + F3; and "其他" (Other) with a checked checkbox for "显示悬浮小球" (Show Floating Ball) and a "悬浮球透明度" (Floating Ball Transparency) slider set to 0%. A "恢复默认" (Reset Default) button is at the bottom left of the settings dialog, and a "说明" (Note) at the bottom right states "所有设置将恢复为程序默认参数" (All settings will be restored to the program's default parameters).

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